



To Reorder Call (281) 647-9700
or 1-800-295-BETA

Product Code: 5250

MEDI-WIPE

Surface Disinfectant

- VIRUCIDAL* • BACTERICIDAL • PSEUDOMONACIDAL • FUNGICIDAL •
- DISINFECTS, DEODORIZES, CLEANS •

GREAT FOR:

- HOSPITALS • MEDICAL OFFICES • DENTAL OFFICES • VETERINARY FACILITIES •
- JANITORIAL, KITCHEN AND BATHROOM FACILITIES •

DIRECTIONS

It is a violation of U.S. Federal Law to use this product in a manner inconsistent with its labeling. This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection. This product is not for use on the inhalation therapy and/or pulmonary diagnostic equipment.

DISINFECTION: Thoroughly wet pre-cleaned, hard, non-porous surface with a wipe, keep wet for 2 minutes (5 minutes if fungus is suspected), and allow to air dry. Use as many wipes as needed for the treated surface to remain wet for the entire contact time. For use as a general disinfectant wipe on hard, non-porous, inanimate surfaces such as chairs, baskets, counter tops, exam tables, cages, toilets, basins, tubs, bed springs, frames, bedside and office furniture, chart racks and equipment, cabinets, door knobs, office machinery, sinks, mirrors, light fixtures and switches, supply and trash carts, telephones, waste cans and laundry hampers, wheel chairs, and walkers, urinals, plastic and other similar hard surfaces. You may use a wipe to pre-clean the surface to be disinfected, but thorough re-wiping with a fresh wipe is required for disinfection. Use on food contact surfaces must be followed by a potable water rinse.

This product is effective against the the following pathogens in 2 minutes at 20°C/68°F, on hard, non-porous surfaces when used as directed: *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Salmonella enterica*, *Escherichia coli* 0157:H7, *Herpes simplex virus 2* (genital Herpes virus), *HIV-1* (AIDS virus)**(At an exposure time of 1 min.), *Mycobacterium tuberculosis*, Methicillin Resistant *Staphylococcus aureus* (MRSA), Vancomycin Resistant *Enterococcus faecalis* (VRE), *Human Hepatitis B virus*, *Human Hepatitis C virus*, *Influenza A* (A2 Japan).

**See back panel for complete HIV-1 (AIDS virus) information.

This product is effective against the the following pathogens in 5 minutes at 20°C/68°F, on hard, non-porous surfaces when used as directed: *Trichophyton mentagrophytes* (Athlete's foot fungus).

DEODORIZING: Wipe hard, non-porous surface thoroughly to eliminate odors from bathrooms, fecal matter, urine, burns, smoke, mildew and mold, decay, stale air, tobacco, sickrooms and musty closets.

CLEANING: To clean surfaces, wipe thoroughly with a fresh wipe, which may be used to deodorize at the same time as cleaning. When used to clean a hard non-porous surface prior to disinfection, a second wipe must be used for disinfection.

THE ACTIVITY OF MEDI-WIPE SURFACE DISINFECTANT AGAINST THE HIV-1 (AIDS VIRUS):

MEDI-WIPE SURFACE DISINFECTANT completely inactivates the HIV-1 (AIDS virus) on hard, non-porous surfaces when used at a contact time of 1 minute at room temperature (20°-25°C). MEDI-WIPE SURFACE DISINFECTANT kills HIV-1 on pre-cleaned environmental surfaces/objects soiled with blood/body fluids, in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type-1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of MEDI-WIPE SURFACE DISINFECTANT. **Contact Time:** HIV-1 virus is inactivated in 1 minute, however, other pathogens require 5 minutes exposure for complete kill. **See disinfection directions. Disposal of Infectious Materials:** Materials that may contain blood or other body fluids should be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal.

ACTIVE INGREDIENTS:

N-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) Dimethyl Benzyl Ammonium Chloride. 0.12%
N-Alkyl (68% C12, 32% C14) Dimethyl Ethyl Benzyl Ammonium Chloride. 0.12%
Isopropyl Alcohol. 41.58%

OTHER INGREDIENTS+: 58.18%

TOTAL: 100.00%

+ Does not include weight of dry wipe.

BETA DATA

FORM:	Premoistened Towel
ODOR:	Lemon/Alcohol
COLOR:	White Cloth
DETERGENCY:	N/A
TOXICITY:	N/A
WETTING ABILITY:	Excellent
PHOSPHATES:	N/A
STORAGE STABILITY:	1 year +
COLD STABILITY:	32°F

FLASH POINT:	73°F TCC
BOILING POINT:	N/A
SOLUBILITY IN WATER:	100%
SPECIFIC GRAVITY:	0.91
% VOLATILE BY VOLUME:	N/A
EVAPORATION RATE:	N/A
BIODEGRADABLE:	Yes
FLAME EXTENSION:	N/A
pH:	5.25

BETA TECHNOLOGY, INC.
CORPORATE HEADQUARTERS
HOUSTON, TX 77218